

Separation of dioxins




Figure: GC/MSD/EISIM TIC showing the isolation of 2,3,7,8-substituted-, tetra-, penta-, and hexa-chlorodibenzo-p-dioxins. One ion (M+2 or M+4) for each congener was monitored (321.9 for TCDD, 355.9 for P5CDD, and 391.9 for H6CDD).

Column: LC-50, 10m x 0.10mm x 0.10 μ m film thickness. Part No. 6010 10010.

Chromatographic Conditions: Oven Temperature 100 $^{\circ}$ C, programmed to 170 $^{\circ}$ C @ 40 $^{\circ}$ C/min, then to 270 $^{\circ}$ C @ 3 $^{\circ}$ C/min, hold 10 minutes at 270 $^{\circ}$ C.

Upper: 2,3,7,8-TCDD at 9.14 minutes (peak 1)

Middle: 1,2,3,7,8-P5CDD at 12.17 minutes (peak 2)

Lower: 1,2,3,4,7,8-H6CDD and 1,2,3,6,7,8-H6CDD at 15.55 minutes (peak 3) and 1,2,3,7,8,9-H6CDD at 17.76 minutes (peak 4).